

Test Name

ACUTE HEPATITIS PANEL

Synonyms**Test Code**

683

Specimen Type

Serum

Preferred Container

Serum separator tube (SST)

Alternate Container

red

Minimum Volume**Temperature/Stability**

Room Temperature: 24 hours or Refrigerated: 72 hours

Handling

Must be centrifuged within 2 hours of collection

Reject

Plasma, Grossly Hemolyzed

Panel Includes

HEP A IgM, HEP B CORE IgM, HBSAG, HCV w/Reflex. If HCV is positive, PCR is performed by Pacific Diagnostic Laboratories.

Rev. 05Sep2024