

**Test Name**

BILIRUBIN, TOTAL AND DIRECT

**Synonyms**

BILIRUBIN PANEL

**Test Code**

112

**Specimen Type**

Serum or plasma

**Preferred Container**

Serum separator tube (SST)

**Alternate Container**

Red or green (heparin)

**Minimum Volume****Temperature/Stability**

Refrigerated: 3 days

**Handling**

Protect from light

**Reject**

Hemolysis

**Panel Includes**

TOTAL BILIRUBIN, DIRECT BILIRUBIN