

Test Name

CREATININE CLEARANCE

Synonyms**Test Code**

987

Specimen Type

24-hour urine **and** serum/plasma

Preferred Container

Urine: >10 mL of entire 24-hour urine collection in clean, leakproof container. Serum: SST

Alternate Container

Urine: 6N HCL or boric acid are acceptable. Serum: green (heparin)

Minimum Volume**Temperature/Stability**

Refrigerated: 6 days

Handling

Provide Total Volume of 24-hour urine. Collect the serum at the time of urine collection or when urine is brought to draw station.

Reject

Grossly hemolyzed serum/plasma. Serum collected not within 24 hours of urine collection.

Panel Includes

CREATININE, EGFR, URINE CREATININE, CREATININE CLEARANCE CALCULATION