

Test Name

LH+FSH+PRL+DHEAS+TSH+TESTO

Synonyms**Test Code**

852

Specimen Type

Serum, morning fasting

Collect 3-4 hours after awakening for PRL

Preferred Container

Serum separator tube (SST)

Alternate Container

red

Minimum Volume**Temperature/Stability**

Room temperature: 8 hours or Refrigerated: 48 hours

Handling

Transfer serum to transport vial.

Reject

Grossly hemolyzed.

Panel Includes

LH+FSH+PRL+DHEAS+TSH+TESTO.