

Test Name

TESTOSTERONE PANEL

Synonyms

FREE AND TOTAL TESTOSTERONE

Test Code

856

Specimen Type

Serum

Preferred: morning specimen

Preferred Container

Serum separator tube (SST)

Alternate Container

Red or green (heparin)

Minimum Volume**Temperature/Stability**

Room temperature: 8 hours or Refrigerated: 48 hours

Handling

Transfer serum to transport vial

Reject

EDTA plasma

Panel IncludesTOTAL TESTOSTERONE, SHBG; FREE AND BIOAVAILABLE TESTOSTERONE
CALCULATIONS