

Test Name

TSH W/ REFLEX TO FREE T4

Synonyms

TSH W/ REFLEX FT4

Test Code

859

Specimen Type

Serum or plasma

For 72 hours before collection, do not take supplements containing biotin

Preferred Container

Serum separator tube (SST)

Alternate Container

red or green (heparin)

Minimum Volume**Temperature/Stability**

Room temperature: 8 hours or Refrigerated: 7 days

Handling**Reject**

Plasma other than heparin

Panel Includes

TSH AND REFLEX TO FREE T4 IF TSH IS ABNORMAL