

Test Name

CALCIUM, CORRECTED (ADJUSTED)

Synonyms**Test Code**

187

Specimen Type

Serum or plasma

Preferred Container

Serum separator tube (SST)

Alternate Container

red or green (heparin)

Minimum Volume**Temperature/Stability**

Room temperature: 7 days or Refrigerated: 22 days

Handling

Separate serum or plasma from cells as soon as possible

Reject

Anticoagulants other than heparin; grossly hemolyzed

Panel Includes

TOTAL CALCIUM, ALBUMIN, CORRECTED CALCIUM CALCULATION