

**Test Name**

CT/GC/TRICH, THINPREP

**Intended Use / Indications**

Detect *Chlamydia trachomatis* and/or *Neisseria gonorrhoeae* and/or *Trichomonas vaginalis*

**Methodology**

Nucleic acid amplification test (NAAT) for the detection of rRNA

**Test Code**

1178

**Specimen Type**

ThinPrep liquid Pap specimen

**Preferred Container**

Entire ThinPrep vial, secured tightly

**Temperature/Stability**

Store 2-30°C: 30 days

**Reject**

SurePath, other specimens

**Panel Includes**

CHLAMYDIA TRACHOMATIS, NEISSERIA GONORRHOEAE, TRICHOMONAS VAGINALIS

**Limitations and Additional Information**

The effects of tampon use or douching have not been assessed. Performance characteristics are available upon request.

Expected TAT: 3 days