

**Test Name**

GLOBULIN WITH A/G RATIO

**Synonyms****Test Code**

149

**Specimen Type**

Serum or plasma

**Preferred Container**

Serum separator tube (SST)

**Alternate Container**

red or green (heparin)

**Minimum Volume****Temperature/Stability**

Room temperature: 1 week or Refrigerated: 1 month

**Handling****Reject**

Plasma other than heparin; grossly hemolyzed

**Panel Includes**

TOTAL PROTEIN, ALBUMIN, GLOBULIN CALCULATION, A/G RATIO