

Test Name

GTT, 2-HOUR, POSTPRANDIAL

Synonyms

2H PP

Test Code

136

Specimen Type

Serum or plasma

Preferred Container

Serum separator tube (SST) or grey (sodium fluoride)

Alternate Container

red or green (heparin)

Minimum Volume**Temperature/Stability**

Room temperature: 8 hours or Refrigerated: 72 hours. Fluoride tube is stable for up to 72 hours at room temperature

Handling

Separate serum/plasma from cells as soon as possible.

Reject

Grossly hemolyzed; blood stored overnight on clot