

Test Name

GTT, OB, 4 specimens (100 grams)

Synonyms**Test Code**

1017

Specimen Type

Serum or plasma, 4 specimens: fasting, 1-hour, 2-hour, 3-hour

Preferred Container

Serum separator tube (SST) or grey (sodium fluoride)

Alternate Container

red or green (heparin)

Minimum Volume**Temperature/Stability**

Room temperature: 8 hours or Refrigerated: 72 hours. Fluoride tube is stable for up to 72 hours at room temperature

Handling

Separate serum/plasma from cells as soon as possible. Label tubes with collection time interval (i.e., fasting, one-hour, two-hour and three-hour).

Reject

Grossly hemolyzed; blood stored overnight on clot

Panel Includes

GLUCOSE x4 (fasting, 1-hour, 2-hour, 3-hour)