

**Test Name**

GTT, OB (50 grams, 1-hour, 1 specimen)

**Synonyms****Test Code**

1016

**Specimen Type**

Serum or plasma

**Preferred Container**

Serum separator tube (SST) or grey (sodium fluoride)

**Alternate Container**

red or green (heparin)

**Minimum Volume****Temperature/Stability**

Room temperature: 8 hours or Refrigerated: 72 hours. Fluoride tube is stable for up to 72 hours at room temperature

**Handling**

Separate serum/plasma from cells as soon as possible.

**Reject**

Grossly hemolyzed; blood stored overnight on clot