

Test Name

OBSTETRIC PANEL WITH HIV

Synonyms

PRENATAL W/ HIV

Test Code

861

Specimen Type

Whole blood and Serum

Preferred Container

Two lavenders (EDTA) and one SST

Alternate Container

red for serum

Minimum Volume**Temperature/Stability**

Room temperature: 8 hours or Refrigerated: 48 hours

Handling

It is preferable to send unopened SST or red-top tube directly to the laboratory

Reject

See individual tests

Panel Includes

OBSTETRIC PANEL PLUS HIV 1/O/2 ANTIBODY/ANTIGEN w/REFLEX.