

Test Name

TRICHOMONAS VAGINALIS, THINPREP

Synonyms

TRICH, T. vaginalis

Intended Use / Indications

Detect *Trichomonas vaginalis* to aid in the diagnosis of trichomoniasis

Methodology

Nucleic acid amplification test (NAAT) for the detection of rRNA

Test Code

1106

Specimen Type

ThinPrep liquid Pap specimen

Preferred Container

Entire ThinPrep vial, secured tightly

Temperature/Stability

Store 2-30°C; 30 days

Reject

SurePath, other specimens

Limitations and Additional Information

Performance has not been evaluated in women less than 14 years of age. Performance characteristics available upon request.

Expected TAT: 3 days