

Test Name

URINE DRUG SCREEN W/ ETOH

Synonyms

URINE DRUG SCREEN WITH ETHANOL; URINE ABUSED DRUG SCREEN

Test Code

984

Specimen Type

Urine

Preferred Container

Clean plastic container with no preservatives and with minimal dead air space.

Alternate Container**Minimum Volume****Temperature/Stability**

Refrigerated: 3 days

Handling

Keep sample tube tightly closed to prevent evaporation of alcohol.

Reject

Serum, urine sample left open

Panel Includes

AMPHETAMINES 300, BARBITURATES 200, BENZODIAZEPINES 200, COCAINE METABOLITE 150, OPIATES 300, OXYCODONE 300, THC 50, plus ETHANOL (QUANTITATIVE, LIMIT OF DETECTION 10 mg/dL).