

Test Name

PROTEIN AND CREATININE, RANDOM URINE

Synonyms

Protein-Creatinine Ratio

Test Code

137

Specimen Type

Urine, random with no preservatives

Preferred Container

Urine: >10 mL of a random collection

Alternate Container**Minimum Volume**

10 mL

Temperature/Stability

Refrigerated: 48 hours

Handling

Collect in a clean plastic container. Label with patient's full name, time of collection, and date of birth.

Reject

Frozen; serum; extremely turbid or bloody urine

Panel Includes

PROTEIN/CREATININE RATIO