# **Test Name**

PAP TEST, CT/GC AND REFLEX HPV

# Synonyms

CERVICAL SCREEN, 21-25 YEARS

# **Test Code**

1186

# Specimen Type

Thinprep liquid Pap specimen

# **Preferred Container**

Entire ThinPrep vial, secured tightly

#### **Alternate Container**

#### Minimum Volume

Temperature/Stability Store 2-30°C

# Handling

For 21-25 years old

#### Reject

SurePath or other specimen types

# Panel Includes

Pap performed by Pacific Diagnostic Laboratories and Reflex to HPV on ASCUS and above (no HPV genotype unless requested); plus CHLAMYDIA/GONOCOCCUS

Rev. 05Nov2019