Test Name

PAP TEST W/ HPV COTEST

Synonyms

CERVICAL SCREEN, 30-64 YEARS

Test Code

1116

Specimen Type

Thinprep liquid Pap specimen

Preferred Container

Entire Thinprep vial, secured tightly

Alternate Container

Minimum Volume

Temperature/Stability

Store 2-30°C

Handling

For 30-64 years old

Reject

Panel Includes

Pap performed by Pacific Diagnostic Laboratories; HPV mRNA E6/E7 and reflex to HPV Genotype if HPV is positive

Rev. 05Nov2019