## **Test Name**

PAP TEST W/ REFLEX HPV

## Synonyms

**CERVICAL SCREEN, 26-29 YEARS** 

## **Test Code**

1117

# Specimen Type

ThinPrep liquid Pap specimen

# **Preferred Container**

Entire Thinprep vial, secured tightly

**Alternate Container** 

### Minimum Volume

Temperature/Stability Store 2-30°C

# Handling

For 26-29 years old

#### Reject

SurePath or other specimen types

### **Panel Includes**

Pap performed by Pacific Diagnostic Laboratories; Reflex to HPV on ASCUS and above (no genotype, unless requested)

Rev. 05Nov2019